


# Wheelchair Lateral Supports

# Wheelchairs

This new design of lateral support offers a number of advantages. The new slimmer design still retains the strength and reliability required to provide consistent, repeatable positioning within a standard wheelchair. They are easy to fit and allow an excellent range of height and width adjustment.

Select the Lateral Pad Size and the tube length. This will show the Lateral Offset from the wheelchair backrest tubes. This gives the maximum distance that is possible and will vary as the width between the lateral supports on each side of the wheelchair is increased or decreased.

## Size Information



Part No.	Lateral Pad Size	LATERAL OFFSET						
C30/17Q	A (13x7.5)	1	3	4	5	8	10.5	13
C30/17S	B (14x7.5)	0.5	2	3	4	7	9.5	12
C30/17T	C (15x9.5)	0.5	2	3	4	7	9.5	12
C30/17U	D (16x9.5)	0.5	2	3	4	7	9.5	12
C30/17V	E (19.5x9.5)	0.5	2	3	4	7	9.5	12
C30/17W	F (25x9.5)	-	-	0.5	1.5	4.5	7	9.5
C30/17X	G (30x9.5)	-	-	0.5	1.5	4.5	7	9.5

All dimensions are in cm.

## YOUR PRESCRIPTION DETAILS

### LATERAL SUPPORT DETAILS

Lateral Pad Size: \_\_\_\_\_

Lateral Offset: \_\_\_\_\_

Quantity: \_\_\_\_\_ pairs

Tube Clamp Size: 3/4": ☐ 7/8": ☐

Fabric Colour: \_\_\_\_\_

Please attach a copy of this sheet with your order.  
See our Fabric Colour Card for further details of upholstery options.

## Compatibility

The Wheelchair Lateral Supports are suitable for use with most wheelchairs that have backrest tubes whose diameter is either 3/4" or 7/8". They are fitted easily using the split clamps and can be adjusted for width. Height adjustment is provided by moving the lateral support on the wheelchair backrest tubes.

## Additional Support

They can be used to provide a symmetrical trunk position where the support offered by the wheelchair itself is not sufficient.

## Kit Contents

The Wheelchair Lateral Support Kit comprises of the following items:

- ☐ Lateral support pad
- ☐ Pad cover
- ☐ 7/8" split or 3/4" standard tube clamp for easy assembly
- ☐ Compression tube
- ☐ All fixings.

## Other Accessories

Lateral Supports are also available separately.

A swingaway version is also available.

LFT029-7

As a company we are committed to promoting a 24 hour approach to postural management through education, design and manufacture. The design of our equipment and training is based on a thorough understanding of the theoretical and practical requirements of the management of posture for people with moderate to severe motor impairment.



**CE** This product is manufactured to comply with the Medical Devices Directive MDD93/42/EEC.



Active Design Ltd  
68K Wyrley Road  
Birmingham B6 7BN  
Tel: (0121) 326 7506  
Fax: (0121) 327 8807  
Web: [www.activedesign.co.uk](http://www.activedesign.co.uk)